

Cheirón[®] 
...we breathe for you.

EU EACH PUMP ✓
TESTED
SER. No 8888 SINCE 1994



VICTORIA

ICONIC SUCTION UNIT

- **EXTREMELY SILENT RUN**
- MAINTENANCE-FREE GERMAN AGGREGATE
- AGGREGATE BODY IS MADE OF ALUMINIUM ALLOY
- RELIABLE AND VIBRATION-FREE RUN
- AGGREGATE IS COOLED BY EXTRA FAN
- DESIGNED FOR CONTINUOUS RUN
- LOW MAINTENANCE REQUIREMENTS
- **6 YEAR WARRANTY ON THE AGGREGATE**
- INNOVATIVE EUROPEAN DESIGN

www.cheiron.eu

- **EXTREMELY SILENT RUN**
- MAINTENANCE-FREE GERMAN AGGREGATE
- AGGREGATE BODY IS MADE OF ALUMINIUM ALLOY
- RELIABLE AND VIBRATION-FREE RUN
- AGGREGATE IS COOLED BY EXTRA FAN
- DESIGNED FOR CONTINUOUS RUN
- LOW MAINTENANCE REQUIREMENTS
- **6 YEAR WARRANTY ON THE AGGREGATE**
- INNOVATIVE EUROPEAN DESIGN



Overflow protection

Vacuum gauge in kPa, mmHg and Bar, colour coded

Illuminated on/off switch meets EU standard EN ISO 10079-1

Ergonomic vacuum regulator

Casing is made of very sturdy, unbreakable, fire retardant material with UV protection

Change-over valve shift the vacuum from one jar to the other with a single move – optional

Jar holder is equipped with suction circuit holder

The stainless eurorail ensures maximum flexibility. You can attach up to four jars of any size and/or any other accessories needed

MAKES LESS NOISE THAN YOUR PC!

A SPECIAL SYSTEM OF NOISE AND VIBRATION DAMPENING FOR AN EXTREMELY QUIET AND SMOOTH OPERATION WITH AN INCREDIBLY SILENT RUN

Non-marking casters with safety brakes.
Antistatic casters - optional

VICTORIA VERSA

The most popular high performance model of Victoria especially suitable for wards, operating rooms, recovery rooms, ICUs, obstetrics and endoscopic applications, etc.

VICTORIA LIPOS

High performance solution designed particularly for liposuction and operating rooms.

VICTORIA PORTABLE

Quiet and effective suction unit for medical practices, wards and endoscopic applications, as well suitable for fast suction in any situation.

VICTORIA THORAX

A suction unit specially designed for common use as well as for thoracic and abdominal surgery drainage and intensive care. The system can be used for long-term up to 30 days continuous drainage.

Thanks to a unique patented solution Victoria Thorax enables full suction output even during drainage.



VICTORIA PORTABLE WITH EURORAIL

Advantage of the unit is that instead of regular holder, unit is equipped with eurorail 25 x 10 mm this allows you to use your holders and jars.



TECHNICAL SPECIFICATIONS

TYPE	PORTABLE	VERSA	THORAX	LIPOS
Suction jar 2 l with press on cover	•	•		•
Suction jar 2 l with screw cover			• 2 pcs	
Suction jar holder for 2 l jar for the side of the unit	•			
Suction jar holder for 2 l jar for eurorail		•	• 2 pcs	•
PVC suction circuit with stop valve	•	•	•	•
Mobile stand		•	•	•
Limiting water valve			•	
Suction capacity nominal +/- 5 %	43 l/min	43 l/min	43 l/min	50 l/min
Suction capacity effective +/- 5 %	38 l/min	38 l/min	38 l/min	45 l/min
Vacuum	0-90 kPa (0-680 mmHg)	0-90 kPa (0-680 mmHg)	0-90 kPa (0-680 mmHg)	0-93 kPa (0-700 mmHg)
Drainage			0 - 40 cmH ₂ O	
Input Power	110 VA	110 VA	110 VA	110 VA
Input Voltage	230V /50-60 Hz and models for 120V ± 10% / 60Hz are also available			
Dimensions (w x h x l) in millimeters	440 x 300 x 310	480 x 900 x 460	480 x 900 x 460	480 x 900 x 460
Dimensions (w x h x l) in inches	20,1 x 11,8 x 13,4	18,1 x 18,9 x 39	18,1 x 18,9 x 39	18,1 x 18,9 x 39
Weight (without accessories)	12 kg (26,5 lb)	21 kg (46,3 lb)	24 kg (52,9 lb)	21 kg (46,3 lb)

OPTIONAL ACCESSORIES

Vacuum gauge in kPa, mmHg and Bar, colour coded



Bottle against overflow Victoria

Order No: 0-52-0002



Support ring for 2 l jar

Order No: 11-5141



Clamp 25x5

Order No: 11-5122



New eurorail for placing a large number of optional accessories

Order No: 11-5104



Disposable secreta bag MONOKIT

Order No: 114-087-001



0,5 – 4 l polycarbonate and polysulfone suction jars, autoclavable

Order No: Price list



MSF microbiological suction filter (20pcs/box)

Order No: 271-022-001



Easy and safe handling



Change-over valve

Order No: 11-5107



IEC power outlet. It is possible to order unit with IEC power outlet



Air-operated on/off footswitch

Order No: 114-047-000



Shelf

Order No: 11-5116



Storage basket

Order No: 11-5105



Holders allow you to place extra jars



Suction circuit holder

Order No: 11-5110

